Anton Paar

pH2Go

Safety Guide

EN	Safety Instructions (Original)	4
DE	Sicherheitshinweise (Übersetzung)	9
FR	Consignes de sécurité (traduction)	14
IT	Istruzioni di sicurezza (traduzione)	19
ES	Instrucciones de seguridad (Traducción)	24
PT	Instruções de segurança (tradução)	29
SV	Säkerhetsanvisningar (översättning)	34
FI	Turvallisuusohjeet (käännös)	39
DA	Sikkerhedsinstruktioner (Oversættelse)	44
NL	Veiligheidsinstructies (vertaling)	49
PL	Instrukcje bezpieczeństwa (tłumaczenie)	54
CS	Bezpečnostní pokyny (překlad)	59
SK	Bezpečnostné pokyny (preklad)	64
HU	Biztonsági útmutató (fordítás)	69
HR	Sigurnosne upute (prijevod)	74
SL	Varnostna navodila (prevod)	79
BG	Инструкции за безопасност (Превод)	84
RO	Instrucțiuni de siguranță (Traducere)	89
LT	Saugos instrukcijos (vertimas)	94
ET	Ohutusjuhend (tõlge)	99
LV	Drošības instrukcijas (tulkojums)	104
EL	Οδηγίες ασφαλείας (Μετάφραση)	109
ZH	安全说明 (译文)	114

SAFETY INSTRUCTIONS (Original)

Anton Paar hereby declares that the radio equipment type pH2Go meets the Directive 2014/53/EU. The full text of the EU Declaration of Conformity is available at the following Internet address: <https://www.anton-paar.com/corp-en/services-support/document-finder/>. Read the safety instructions before using the device. Make this safety instructions easily accessible to all persons working with the device. Follow all tips and instructions in this safety instructions to ensure the correct use and safe functioning of the device.

CONVENTIONS FOR SAFETY MESSAGES

The following conventions for safety messages are used:

WARNING

Warning indicates a hazardous situation, which, if not avoided, could result in death or serious injury.

CAUTION

Caution indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

NOTICE

Notice indicates a situation which, if not avoided, could result in damage to the property.

LIABILITY

- > This safety instruction does not claim to address all safety issues associated with the use of the device and samples. It is your responsibility to establish health and safety practices and to determine the applicability of regulatory limitations.
- > Anton Paar ConsumerTec GmbH warrants the proper functioning of the device only if no modifications have been made to mechanics, electronics, firmware or software.
- > Only use the device for the purpose described in the documentation. Anton Paar ConsumerTec GmbH is not liable for damages caused by the incorrect use of the device.
- > The results delivered by the device depend not only on the correct functioning of the device, but also on various other factors. We therefore recommend you have the results checked (e.g. plausibility tested) by skilled personnel before consequential actions are taken based on the results.

OPERATION IN AREAS WITH RISK OF EXPLOSION

WARNING

- > The device is not explosion-proof and therefore must not be operated in areas with risk of explosion.

INSTALLATION AND USE

WARNING

- > Do not use the device in medical areas.
- > Use the device as a stand-alone unit combined with a mobile phone, do not combine the device with other devices for automation purposes.
- > Do not use any spare parts other than those supplied or approved by Anton Paar ConsumerTec GmbH.
- > In case of damage or malfunction, do not continue operating the device. Do not operate the device under conditions which could result in damage to goods or injuries or loss of life.
- > Make sure all operators are trained to use the device safely and correctly before starting any applicable operations.
- > Ensure the device is sufficiently supervised during operation.
- > Keep the device away from children.

PERSONAL SAFETY

WARNING

- > Always keep a distance of at least 15 cm (6 inch) between the device and an implanted cardiac pacemaker or defibrillator.

CAUTION

- > During operation of the device, the operator must be equipped with appropriate personal protective equipment (safety goggles, gloves).

GENERAL PRECAUTIONS

WARNING

- > Observe and adhere to your national safety regulations regarding the handling of all substances associated with your measurements (e.g. use safety goggles, gloves, respiratory protection, etc.).
- > This device is designed for a wide range of applications in various areas. The use of this device requires experience and knowledge in working with the used substances.
- > Before using substances, always check the corresponding safety data sheet.
- > Fill only substances which will not pose a risk if they leak from the device in case of damage or malfunction.
- > Take measures that spilled liquids cannot get into plug connections or venting slots of electrical appliances.
- > Take care that the liquids (samples and cleaning agents) that you use are chemically compatible when they come into contact with each other. They must not react exothermically or produce harmful substances.
- > Do not leave the device in a location exposed to direct sunlight or near a heat source for an extended period of time.
- > Do not change the ambient temperature of the device abruptly.

**CAUTION**

- > Be careful with hot samples. Handling samples with temperatures more than 50 °C (122 °F) bears the risk of burns. Wear appropriate protective clothing or take other protective measures when you handle hot samples.

NOTICE

- > Before a measurement, check the wetted parts of the instrument for chemical resistance to the samples and cleaning agents used. Do not use samples which will attack the wetted parts.

Precautions for highly flammable samples and cleaning agents**WARNING**

- > Keep potential sources of ignition, like sparks or open flames, at a safe distance from the device.
- > Store only the minimum required amount of sample, cleaning agents and other flammable materials near the device.
- > Do not spill sample/cleaning agents or leave their containers uncovered. Immediately remove spilled sample/cleaning agents.
- > Make sure that the setup location is sufficiently ventilated. The environment of the instrument must be kept free of flammable gases and vapors.
- > Provide a fire extinguisher.

PRECAUTIONS FOR BUFFER SOLUTIONS**CAUTION**

- > It is recommended to read the manufacturer's safety data sheets before first usage.
- > Do not eat, drink, or smoke while handling the buffers. Wash hands thoroughly after use. Do not swallow any buffer liquids.
- > Do not inhale vapors, fumes, or mist. If ventilation is insufficient, use appropriate respiratory protection in accordance with local regulations. Ensure that the work area is well-ventilated.
- > Handle with care to avoid spills and splashes. In case of a spill, absorb the liquid with paper and clean the area with plenty of water.
- > For storage conditions of the buffer solutions it is recommended to read the manufacturer's safety data sheet.
- > Dispose of waste according to local environmental regulations and the manufacturer's safety data sheet.

PRECAUTIONS FOR SENSOR AND DIAPHRAGM**NOTICE**

- > The sensor is designed for pH measurements. Do not use the sensor for other applications as this may cause damage.
- > Do not use the sensor and diaphragm outside the specified temperature range (not higher than 80°C / 180°F) as this might result in failure or damage to the sensor and diaphragm.

- > Samples must be aqueous solutions and compatible with the sensor and diaphragm.
- > Avoid exposure to sticky, oily, or particulate-rich samples that may clog the sensor and diaphragm. Rinse the sensor and diaphragm with distilled water immediately after each use. If clogged, soak the sensor and diaphragm in mild cleaning solution.
- > Never scrub or use abrasive materials on the sensor or diaphragm; this may permanently damage the sensor surface.
- > Do not use hydrofluoric acid, acetone, MEK or similar agents to clean the sensor.
- > The sensor should be calibrated regularly using fresh, clean buffer solutions. Ensure that the buffer solutions are not reused or contaminated.
- > Do not store the sensor and diaphragm in distilled or deionized water.
- > After a long-term storage always revitalize the sensor with a saturated KCl solution.
- > In case of damage or malfunction of sensor, or diaphragm, do not continue operating the device.

BATTERY HANDLING**NOTICE**

- > Turn the device off prior to exchanging the batteries.
- > Make sure the sample area is empty and dry before exchanging the batteries to avoid spillage of liquids.
- > Only open the device in a clean and dry environment.
- > Do not open the compartment containing the electrical components of the device.
- > The inside of the device does not need any cleaning. Do not use liquid cleaning agents on interior parts of the device.
- > Use 1.5 V batteries type AAA.
- > In case you are not using the device for a long time, remove the batteries as they could leak or be damaged.
- > Do not submerge the device in water or any other cleaning agents or other liquids.
- > After replacing the batteries, be sure to only use the original screws that include an O-ring to fix the battery compartment cover. Ensure that the screws are tightened properly so that no liquid can get inside. Do not overtighten the screws. Using excessive force may permanently damage your device.
- > When replacing the batteries, check if the O-ring between the bottom and the top cover is not damaged and in place.
- > Insert only equally charged batteries of the same type.
- > Consider the polarity of the batteries when inserting them into the battery compartment.

**WARNING**

- > Do not leave batteries lying unattended.
- > Keep batteries away from children. They can be dangerous for children if used as toys.
- > If a battery has been swallowed, seek medical advice immediately.
- > Leaking or damaged batteries can cause burns if they come into contact with skin. Use gloves for handling.
- > Never short-circuit open batteries.
- > Do not expose batteries to heat or throw them into fire.
- > Do not charge non-rechargeable batteries, there is a risk of explosion.
- > Do not insert damaged batteries in the battery compartment, they can cause damage to the device and fire.

REPAIRS

NOTICE

- > Repair procedures may be carried out only by authorized persons or by Anton Paar ConsumerTec GmbH.

DEVICE RETURNS

- > For returns, contact the vendor of the device.
- > The device must be cleaned before return.

WARNING

- > You must not return devices which are contaminated by radioactive materials, infectious agents, or other harmful substances that can pose health hazards.

DISPOSAL

- > Concerning the disposal of the device, observe the legal requirements in your country.

Address of the instrument producer:
Anton Paar ConsumerTec GmbH
 Anton-Paar-Str. 20
 A-8054 Graz/Austria – Europe

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