



Certificate of Compliance

EasyDens

Part No 234929

Dear customer,

With this certificate the manufacturer declares that this instrument/ equipment was inspected, also tested, and accepted in accordance with our quality management system. It complies with the specifications according to the above listed part number. The following tests have been successfully performed:

- Optical check of electronic board, surface finish, wiring system.
- Functional test of the instrument, including measuring cell, temperature measurement.
- Test of the interface.
- Measurements: Adjustment with air at 20 °C and distilled water at 10 °C also 20 °C, and 40 °C.
- Check with distilled water and Ethanol at different temperatures.
- Check of all labels.
- Configuration check.
- Final check.